

## Product Specification of Latex Powder Free 14.0 Mil Silverlined

**Name of Product** : Latex Powder Free Diamond Silverlined Gloves.  
**Name of Manufacturer** : Innovative Gloves Co. Ltd. Hatyai, Thailand.  
**Base Material** : Natural Rubber Latex.  
**Characteristics** : Non Sterile, Ambidextrous Beaded Cuff.  
**Available Sizes** : Small (7), Medium (8), Large (9), X-Large (10), XX- Large (11).  
**Surface** : Diamond Textured.  
**Surface Treatment** : Off Line Chlorinated.  
**Colour** : Blue and White.  
**Length** : 300 ± 10 mm.

### Physical Specifications:

Size	Thickness - Single Wall (mm)			Width (mm)	
	Tip		Palm		Cuff
	Smooth Surface	Textured Surface	Textured Surface		Smooth Surface
Small (7)	0.37 ± 0.01	0.66 ± 0.01	0.74 ± 0.01	0.21 ± 0.01	85 ± 5
Medium (8)	0.37 ± 0.01	0.66 ± 0.01	0.74 ± 0.01	0.21 ± 0.01	95 ± 5
Large (9)	0.37 ± 0.01	0.66 ± 0.01	0.74 ± 0.01	0.21 ± 0.01	105 ± 5
X-Large (10)	0.37 ± 0.01	0.66 ± 0.01	0.74 ± 0.01	0.21 ± 0.01	115 ± 5
XX-Large (11)	0.37 ± 0.01	0.66 ± 0.01	0.74 ± 0.01	0.21 ± 0.01	120 ± 3

### Physical Properties:

Physical Properties	Force at Break (N)	Tensile Strength (MPa)	Ultimate Elongation (%)
Before Aging	6 N Min.	21 MPa Min.	700% Min.
After Aging	6 N Min.	18 MPa Min.	500% Min

**Total Residue Content** : < 2.0 milligrams Per Glove.

**Total Extractable Protein Content** : < 50 µg/g Of Glove.

**Packaging (Normal)** : 50 ± 1 pieces, 100 ± 2 pieces In One Dispenser Box,  
10 Dispenser Box In One Carton.

**Labeling On Dispenser (Normal)** : Product Name, Size, Lot No. Quantity,  
Manufacturer's Address, Warning About Latex Allergy.

**Labeling On Carton (Normal)** : Product Name, Size, Lot No. Quantity, Country Of Origin.

