

Product Specification of Nitrile 5.0 mil

Name of Product : Nitrile Powder Free Fully Textured Gloves.
Name of Manufacturer : Innovative Gloves Co. Ltd. Hatyai, Thailand.
Base Material : Nitriles.
Characteristics : Non Sterile, Ambidextrous Beaded Cuff.
Available Sizes : X-Small, Small, Medium, Large, X-Large, XX-Large.
Surface : Fully Textured.
Surface Treatment : Polymer Coating.
Colour : Blue, Black, Orange, Light Blue.
Length : 240 ± 5 mm.

Physical Specifications:

Size	Thickness - Single Wall (mm)			Width (mm)
	Tip	Palm	Cuff	
X-Small	0.14 ± 0.01	0.11 ± 0.01	0.09 ± 0.01	75 ± 5
Small	0.14 ± 0.01	0.11 ± 0.01	0.09 ± 0.01	85 ± 5
Medium	0.14 ± 0.01	0.11 ± 0.01	0.09 ± 0.01	95 ± 5
Large	0.14 ± 0.01	0.11 ± 0.01	0.09 ± 0.01	105 ± 5
X-Large	0.14 ± 0.01	0.11 ± 0.01	0.09 ± 0.01	115 ± 5
XX-Large	0.14 ± 0.01	0.11 ± 0.01	0.09 ± 0.01	120 ± 3

Physical Properties:

Physical Properties	Force at Break (N)	Tensile Strength (MPa)	Ultimate Elongation (%)
Before Aging	6 N Min.	22 MPa Min.	500% Min.
After Aging	6 N Min.	22 MPa Min.	400% Min.

Total Residue Content : < 2.0 milligrams Per Glove.
Packaging (Normal) : 100 ± 2 Pieces In One Dispenser Box,
 10 Dispenser Box In One Carton.
Labeling On Dispenser (Normal) : Product Name, Size, Lot No. Quantity,
 Manufacturer's Address, Warning About Latex Allergy
Labeling On Carton (Normal) : Product Name, Size, Lot No. Quantity, Country Of Origin.

