

## Product Specification of Nitrile 6.0 mil

<b>Name of Product</b>	: Nitrile Powder Free Fully Textured Gloves.
<b>Name of Manufacturer</b>	: Innovative Gloves Co. Ltd. Hatyai, Thailand.
<b>Base Material</b>	: Nitriles.
<b>Characteristics</b>	: Non Sterile, Ambidextrous Beaded Cuff.
<b>Available Sizes</b>	: X-Small, Small, Medium, Large, X-Large, 2X-Large, 3X-Large.
<b>Surface</b>	: Fully Textured.
<b>Surface Treatment</b>	: Polymer Coating.
<b>Colour</b>	: Black, Blue, Orange, Royal Purple.
<b>Length</b>	: 240 ± 5 mm.

### Physical Specifications:

Size	Thickness - Single Wall (mm)			Width (mm)
	Tip	Palm	Cuff	
X-Small	0.16 ± 0.01	0.13 ± 0.01	0.10 ± 0.01	75 ± 5
Small	0.16 ± 0.01	0.13 ± 0.01	0.10 ± 0.01	85 ± 5
Medium	0.16 ± 0.01	0.13 ± 0.01	0.10 ± 0.01	95 ± 5
Large	0.16 ± 0.01	0.13 ± 0.01	0.10 ± 0.01	105 ± 5
X-Large	0.16 ± 0.01	0.13 ± 0.01	0.10 ± 0.01	115 ± 5
2X-Large	0.16 ± 0.01	0.13 ± 0.01	0.10 ± 0.01	120 ± 3
3X-Large	0.16 ± 0.01	0.13 ± 0.01	0.10 ± 0.01	127 ± 3

### Physical Properties:

Physical Properties	Force at Break (N)	Tensile Strength (MPa)	Ultimate Elongation (%)
Before Aging	6 N Min.	22 MPa Min.	500% Min.
After Aging	6 N Min.	22 MPa Min.	400% Min.

**Total Residue Content** : < 2.0 milligrams Per Glove.

**Packaging (Normal)** : 100 ± 2 pieces In One Dispenser Box,  
10 Dispenser Box In One Carton.

**Labeling On Dispenser (Normal)** : Product Name, Size, Lot No. Quantity,  
Manufacturer's Address, Warning About Latex Allergy.

**Labeling On Carton (Normal)** : Product Name, Size, Lot No. Quantity, Country Of Origin.

