

Product Specification of Latex Powder Free Gloves

Name of Product : Latex Powder Free Fully Textured Polymer Coated Gloves.
Name of Manufacturer : Innovative Gloves Co. Ltd. Hatyai, Thailand.
Base Material : Natural Rubber Latex.
Characteristics : Non Sterile, Ambidextrous Beaded Cuff.
Available Sizes : X -Small, Small, Medium, Large, X-Large, XX-Large.
Surface : Fully Textured.
Surface treatment : Off Line Chlorinated.
Colour : White.
Length : 240 ± 5 mm.

Physical Specifications:

Size	Thickness - Single Wall (mm)			Width (mm)
	Tip	Palm	Cuff	
X-Small	0.14 ± 0.01	0.11 ± 0.01	0.09 ± 0.01	75 ± 5
Small	0.14 ± 0.01	0.11 ± 0.01	0.09 ± 0.01	85 ± 5
Medium	0.14 ± 0.01	0.11 ± 0.01	0.09 ± 0.01	95 ± 5
Large	0.14 ± 0.01	0.11 ± 0.01	0.09 ± 0.01	105 ± 5
X-Large	0.14 ± 0.01	0.11 ± 0.01	0.09 ± 0.01	115 ± 5
XX-Large	0.14 ± 0.01	0.11 ± 0.01	0.09 ± 0.01	120 ± 3

Physical Properties:

Physical Properties	Force at Break (N)	Tensile Strength (MPa)	Ultimate Elongation (%)
Before Aging	6 N Min.	21 MPa Min.	700% Min.
After Aging	6 N Min.	18 MPa Min.	500% Min

Total Residue Content : < 2.0 milligrams Per Glove.

Total Extractable Protein Content : < 50 µg/g Of Glove.

Packaging (Normal) : 100 ± 2 Pieces In One Dispenser Box,
 10 Dispenser Box In One Carton.

Labeling On Dispenser (Normal) : Product Name, Size, Lot No. Quantity,
 Manufacturer's Address, Warning About Latex Allergy.

Labeling On Carton (Normal) : Product Name, Size, Lot No. Quantity, Country Of Origin.

