

Product Specification of Nitrile 9.0 Mil Micro Diamond Textured

Name of Product : Nitrile Powder Free Micro Diamond Textured Gloves.
Name of Manufacturer : Innovative Gloves Co. Ltd. Hatyai, Thailand.
Base Material : Nitriles.
Characteristics : Non Sterile, Ambidextrous Beaded Cuff.
Available Sizes : Small, Medium, Large, X-Large, 2X- Large, 3X-Large.
Surface : Micro Diamond Textured.
Surface Treatment : Polymer Coating.
Colour : Royal Blue.
Length : 290 ± 5 mm.

Physical Specifications:

Size	Thickness - Single Wall (mm)				Width (mm)
	Tip		Palm	Cuff	
	Smooth Surface	Textured Surface	Textured Surface	Smooth Surface	
Small	0.14 ± 0.01	0.34 ± 0.02	0.30 ± 0.02	0.14 ± 0.01	85 ± 5
Medium	0.14 ± 0.01	0.34 ± 0.02	0.30 ± 0.02	0.14 ± 0.01	95 ± 5
Large	0.14 ± 0.01	0.34 ± 0.02	0.30 ± 0.02	0.14 ± 0.01	105 ± 5
X Large	0.14 ± 0.01	0.34 ± 0.02	0.30 ± 0.02	0.14 ± 0.01	115 ± 5
2X Large	0.14 ± 0.01	0.34 ± 0.02	0.30 ± 0.02	0.14 ± 0.01	120 ± 3
3X Large	0.14 ± 0.01	0.34 ± 0.02	0.30 ± 0.02	0.14 ± 0.01	127 ± 3

Physical Properties:

Physical Properties	Force at Break (N)	Tensile Strength (MPa)	Ultimate Elongation (%)
Before Aging	6 N Min.	22 MPa Min.	500% Min.
After Aging	6 N Min.	22 MPa Min.	400% Min.

Total Residue Content : < 2.0 milligrams Per Glove.

Packaging (Normal) : 50 ± 1 or 100 ± 2 pieces In One Dispenser Box,
 10 Dispenser Box In One Carton.

Labeling On Dispenser (Normal) : Product Name, Size, Lot No. Quantity,
 Manufacturer's Address, Warning About Latex Allergy.

Labeling On Carton (Normal) : Product Name, Size, Lot No. Quantity, Country Of Origin.



CE0321

