

Product Specification of Nitrile 9.0 Mil Diamond Textured Gloves

Name of Product : Nitrile Powder Free Heavy Duty Gloves.
Name of Manufacturer : Innovative Gloves Co. Ltd. Hatyai, Thailand.
Base Material : Nitriles.
Characteristics : Non Sterile, Ambidextrous Beaded Cuff.
Available Sizes : Small(7), Medium(8), Large(9), X-Large(10), XX-Large(11).
Surface Treatment : Off Line Chlorinated.
Surface : Diamond Textured.
Colour : Dark Blue.
Length : 300 ± 10 mm.

Physical Specifications:

Size	Thickness - Single Wall (mm)				Width (mm)
	Tip		Palm	Cuff	
	Smooth Surface	Textured Surface	Textured Surface	Smooth Surface	
Small(7)	0.24 ± 0.01	0.40 ± 0.02	0.45 ± 0.02	0.12 ± 0.01	85 ± 5
Medium(8)	0.24 ± 0.01	0.40 ± 0.02	0.45 ± 0.02	0.12 ± 0.01	95 ± 5
Large(9)	0.24 ± 0.01	0.40 ± 0.02	0.45 ± 0.02	0.12 ± 0.01	105 ± 5
X Large(10)	0.24 ± 0.01	0.40 ± 0.02	0.45 ± 0.02	0.12 ± 0.01	115 ± 5
XX Large(11)	0.24 ± 0.01	0.40 ± 0.02	0.45 ± 0.02	0.12 ± 0.01	120 ± 3

Physical Properties:

Physical Properties	Force at Break (N)	Tensile Strength (MPa)	Ultimate Elongation (%)
Before Aging	6 N Min.	22 MPa Min.	500% Min.
After Aging	6 N Min.	22 MPa Min.	400% Min.

Total Residue Content : < 2.0 milligrams Per Glove.

Packaging (Normal) : 50 ± 1 pieces In One Dispenser Box,
10 Dispenser Box In One Carton.

Labeling On Dispenser (Normal) : Product Name, Size, Lot No. Quantity,
Manufacturer's Address, Warning About Latex Allergy.

Labeling On Carton (Normal) : Product Name, Size, Lot No. Quantity, Country Of Origin.

